


Case Report

Prostate MRI

MRI-targeted
TRUS-guided biopsy



Olea Medical
93, avenue des Sorbiers
13600 La Ciotat - FRANCE
TEL + 33 (0)4 42 71 24 20
FAX +33 (0)4 42 71 24 27
contact@olea-medical.com


Olea Medical USA
1955 Massachusetts Avenue
Suite 12 & 14
02140 Cambridge, MA

Patient history

Sixty-four y/o man presenting with a raising PSA level. The more recent value is 7ng/ml. Digital rectal examination is non-suspicious and prostate volume is estimated at 30ml. A multiparametric MRI is performed prior to TRUS-guided biopsies. The MRI protocol includes three sequences : one T2 weighted sequence (T2W), one diffusion weighted sequence (DW) and one dynamic contrast enhanced sequence (DCE). The T2 weighted sequence is a 3D acquisition with a voxel size of 0.35cc (0.7 x 0.7 x 0.7mm). The DW sequence is acquired with several b-values (b100-200-400-1000) from which is extracted the Apparent Diffusion Coefficient value. The DCE sequence is a gradient echo sequence after IV bolus injection (0.1mmol/kg) of a chelate of gadolinium (Gadovist, Bayer). Several phases are acquired with a temporal resolution of 8.5s for a total acquisition time of 5 mn (35 phases).

Post-processing is performed on a dedicated workstation (Olea Sphere™, Olea Medical®, La Ciotat, France) which allows for the measurement of dynamic quantitative parameters, extracted from the Tofts model (1). K_{trans} or transfer constant corresponds to the wash-in. K_{ep} is the rate constant and corresponds to the wash-out. The area under the concentration of gadolinium 60s after the beginning of contrast injection (AUGC), corresponding to the maximum peak enhancement, is also calculated. Absolute values of the three parameters are color-coded and displayed on screen.


Mp-MRI (Fig. 1) shows two suspicious areas within the peripheral zone corresponding to two significant tumors diagnosed by MRI-targeted TRUS-guided biopsies (Fig.2) with TRUS-MRI image fusion (Koealis system, Grenoble, France) and confirmed at pathological examination of the radical prostatectomy specimen (Fig. 3).



■ Figure 1a & 1b: multiparametric MRI.

1a: One lesion is visible in the left apex (arrow). It has a discrete low signal intensity on the T2 image, but impeded diffusion is clearly visible, showing a focal low signal on the ADC map and a high intensity area on the long b-value (B1000). DCE-MRI shows a bright area corresponding to high values of K_{trans} , K_{ep} and AUGC. Note the physiological hypervascularity of the TZ (*).

1b: At the right prostate base, a second lesion is present (arrow), barely visible on T2W-MRI and on the ADC map, but highly conspicuous on trace diffusion images. DCE-MRI shows an increased permeability of tumor vessels on source DCE images (white arrow, perm), on the right postero-lateral aspect of the prostate, but color coded images of quantitative parameters (K_{trans} and AUGC) extracted from the post-processing of the dedicated workstation (OleaSphere™) make the diagnosis obvious.




■ Figure 2: MRI targeted TRUS guided biopsies using the Koealis system which allows for a TRUS-MRI image registration after deformation of the MRI prostate shape to accurately adjust it to the TRUS prostate shape. A region of interest (ROI, yellow tag) has been placed on the MRI targets. Once image fusion has been performed, the ROI's can be displayed indifferently on the MRI or the TRUS image, allowing for an accurate targeting of the lesion on the basis of the MRI findings. Three biopsies have been performed in the apical lesion (1-3, green tags) and two in the lesion of the left base (4-5). Bilateral Gleason score 7 (3+4) cancer. L: left side. R: right side

Discussion.

The ability of multiparametric MRI to localise tumor foci of prostate cancer is now established (2). DW- and DCE-MRI have increased the performance of T2W-MRI (3), which lacks specificity. DW-MRI shows the impeded diffusion of areas with high cellular density, which characterizes tumor foci of prostate cancer, especially when they contain high Gleason grades (4). DCE-MRI with dedicated workstations, like the OleaSphere™ system, enables a fast and accurate reading of the different quantitative (K_{trans}, K_{ep}, v_e, v_p) parameters (1) depending on the type of pharmacokinetic model used. Color coded maps make tumor foci highly conspicuous and this capacity can also be used to assess kinetics of Gadolinium with semi-quantitative parameters, like wash-in, wash-out, time to peak or maximum peak enhancement (5). These parameters can be assessed visually (6, 7) on the source images thanks to the different color coded maps of semi-quantitative parameters available on the OleaSphere™ workstation (8). These color coded images can be combined with placement of regions of interest (ROI's) to obtain different types of curve of enhancement (9).

References

- Tofts PS, Brix G, Buckley DL, Evelhoch JL, Henderson E, Knopp MV, et al. Estimating kinetic parameters from dynamic contrast-enhanced T1-weighted MRI of a diffusible tracer: standardized quantities and symbols. *J Magn Reson Imaging*. 1999 Sep;10(3):223-32.
- Cornud F, Delongchamps NB, Mozer P, Beuvon F, Schull A, Muradyan N, et al. Value of multiparametric MRI in the work-up of prostate cancer. *Current urology reports*. [Review]. 2012 Feb;13(1):82-92.
- Delongchamps NB, Rouanne M, Flam T, Beuvon F, Liberatore M, Zerbib M, et al. Multiparametric magnetic resonance imaging for the detection and localization of prostate cancer: combination of T2-weighted, dynamic contrast-enhanced and diffusion-weighted imaging. *BJU international*. [Evaluation Studies]. 2011 May;107(9):1411-8.
- Afaq A, Koh DM, Padhani A, van As N, Sohaib SA. Clinical utility of diffusion-weighted magnetic resonance imaging in prostate cancer. *BJU international*. [Research Support, Non-U.S. Gov't Review]. 2011 Dec;108(11):1716-22.
- Kiessling F, Lichy M, Grobholz R, Heilmann M, Farhan N, Michel MS, et al. Simple models improve the discrimination of focal abnormalities from the peripheral gland by T1-weighted dynamic MRI. *Eur Radiol*. 2004 Oct;14(10):1793-801.
- Villers A, Puech P, Mouton D, Leroy X, Ballereau C, Lemaître L. Dynamic contrast enhanced, pelvic phased array magnetic resonance imaging of localized prostate cancer for predicting tumor volume: correlation with radical prostatectomy findings. *J Urol*. 2006 Dec;176(6 Pt 1):2432-7.
- Rouviere O, Papillard M, Girouin N, Boutier R, Rabilloud M, Riche B, et al. Is it possible to model the risk of malignancy of focal abnormalities found at prostate multiparametric MRI? *European radiology*. 2012 May;22(5):1149-57.
- Frankel T, Hamm B, Hricak H. Dynamic contrast-enhanced magnetic resonance imaging and pharmacokinetic models in prostate cancer. *European radiology*. [Research Support, Non-U.S. Gov't Review]. 2011 Mar;21(3):616-26.



■ Figure 4 : Pathological examination of the radical prostatectomy specimen. The main (index) lesion (1.6cc) is on the left side (arrows, a), at the apex. The contralateral tumor, located in the prostate base, has a significant volume (0.8cc). Bilateral Gleason score 7 lesion with 40% grade.

Our case illustrates how mp-MRI can identify men with potentially clinical bilateral significant PCa prior to biopsy. It also shows that, once a target has been identified on MRI, the physician must have this information available at the time of biopsy to match as accurately as possible the needle tract and the target, hence the concept of TRUS-MRI image registration (or image fusion) to plan and to guide the biopsy. An accurate TRUS-MRI fusion system has to take into account the difference in prostate shape between TRUS and MRI and also consider patient and prostate motion during TRUS examination. The Koelis system (Koelis Inc., La Tronche, France) has this capability and allows, as illustrated in our case, for a precise MRI-targeting of focal abnormalities during TRUS guided biopsies (10).

In conclusion, our case suggests that the use of MRI in prostate cancer management is very promising, not only in patients with a diagnosed PCa, but also for men before prostate biopsy. Standardization of criteria which define a target is mandatory. Different scoring systems are currently used (3, 9, 11, 12). All include an individual three or five-point scale T2/DW/DCE score and an overall five point scale Likert score with a trend to take into account the dominant sequence according to the zone of origin of the lesion (DWI-MRI for the PZ and T2W-MRI for the TZ). Because the degree of suspicion on MRI is a powerful predictor of significant cancer (13-16), a standardized report is crucially required to define what is benign or probably benign (score 1-2), equivocal (score 3) and probably malignant or malignant (score 4-5).

- Barentsz JO, Richenberg J, Clements R, Choyke P, Verma S, Villeirs G, et al. ESUR prostate MR guidelines 2012. *European radiology*. 2012 Apr;22(4):746-57.
- Cornud F, Brois L, Delongchamps NB, Portalez D, Malavaud B, Renard-Penna R, et al. TRUS-MRI image registration: a paradigm shift in the diagnosis of significant prostate cancer. *Abdominal imaging*. 2013 Jul 17.
- Arumainayagam N, Ahmed HU, Moore CM, Freeman A, Allen C, Sohaib SA, et al. Multiparametric MR Imaging for Detection of Clinically Significant Prostate Cancer: A Validation Cohort Study with Transperineal Template Prostate Mapping as the Reference Standard. *Radiology*. 2013 Apr 5.
- Natarajan S, Marks LS, Margolis DJ, Huang J, Macairan ML, Lieu P, et al. Clinical application of a 3D ultrasound-guided prostate biopsy system. *Urol Oncol*. 2011 May-Jun;29(3):334-42.
- Portalez D, Mozer P, Cornud F, Renard-Penna R, Misrai V, Thoulouzan M, et al. Validation of the European Society of Urogenital Radiology Scoring System for Prostate Cancer Diagnosis on Multiparametric Magnetic Resonance Imaging in a Cohort of Repeat Biopsy Patients. *European urology*. 2012 Jun 27.
- Sonn GA, Chang E, Natarajan S, Margolis DJ, Macairan M, Lieu P, et al. Value of Targeted Prostate Biopsy Using Magnetic Resonance-Ultrasound Fusion in Men with Prior Negative Biopsy and Elevated Prostate-specific Antigen. *European urology*. 2013 Mar 17.
- Sonn GA, Natarajan S, Margolis DJ, Macairan M, Lieu P, Huang J, et al. Targeted biopsy in the detection of prostate cancer using an office based magnetic resonance ultrasound fusion device. *The Journal of urology*. [Research Support, N.I.H., Extramural Research Support, Non-U.S. Gov't]. 2013 Jan;189(1):86-91.
- Rud E, Bao E, Eggesbo HB. MRI and Ultrasound-guided Prostate Biopsy Using Soft Image Fusion. *Anticancer research*. 2012 Aug;32(8):3383-9.



Olea Medical
93, avenue des Sorbiers
13600 La Ciotat - FRANCE
TEL + 33 (0)4 42 71 24 20
FAX +33 (0)4 42 71 24 27
contact@olea-medical.com

Olea Medical USA
1955 Massachusetts Avenue
Suite 12 & 14
02140 Cambridge, MA

