



Jan Casselman

Professeur de Radiologie

AZ St. Jan Brugge-Oostende AV, Bruges, Belgique

Selon vous, quelles sont les spécificités essentielles d'une solution de post-traitement pour un radiologue ?

"Once the often complex and numerous ENT MR and/or CT images are acquired, dedicated image processing software is needed to analyze these images and to get access to all information imbedded in these images.

The software should speed-up and automate frequently used image analysis and also should open the more complex image processing to radiologists who are working in a non-research environment. This is exactly what the Olea Sphere® software is offering.

Quels bénéfices la solution de post-traitement Olea Sphere® vous apporte-t-elle dans votre spécialité ?

Unlike other processing software, the Olea Sphere® software can be tailored to the specific needs of every radiologist and different processing steps can be grouped in a "one or two click workflow, making life easy and saving a lot of time.

Registration or matching of different MR sequences of one study (for the evaluation of neurovascular conflicts), same MR sequences made at different follow-up dates (for the evaluation of 7th nerve schwannoma growth) and of MR and CT images (for the exact localization of cholesteatomas in the middle ear and of cochlear implants inside the cochlea) are some of the routine ENT applications which can be performed in a reliable and fast way with the Olea Sphere® software.

However, the Olea Sphere® software becomes even more indispensable when more complex image processing is needed.

**“Unlike other processing software,
the Olea Sphere® software can be tailored
to the specific needs of every radiologist...”**

Multi-parametric advanced ENT tumor imaging today consists of routine anatomical MR images completed by perfusion, permeability and Intra Voxal Incoherent Motion (IVIM) diffusion weighted images, and is crucial in the prediction of treatment response, treatment monitoring and prediction of final outcome.

This multi-parametric approach implies that thousands of images are processed in an acceptable time and that complex calculations as IVIM (resulting in the "true" diffusion, excluding false diffusion caused by vascular flow) are possible.

Moreover, the software again must be able to automate most of the processing steps and must display all the multi-parametric results in a clear easy to interpret custom made layout.

All the above is possible with the Olea Sphere® software and it will even be difficult to find other software which is for instance able to calculate IVIM values.

Therefore the Olea Sphere® software helps you to save time, helps you to get all information out of the images and gives you access to high end image processing which was before only available and affordable in a research environment."